




# HIGH PERFORMANCE 'NEAT-A-WARD'® PATIENT WALLBEDS

FOR MIDWIFE-LED BIRTHING UNITS















# HIGH PERFORMANCE 'NEAT-A-WARD'® PATIENT WALLBEDS











# Instantly Available Rest/Overnight Facilities

Our high performance 'NEAT-A-WARD'® Patient Wallbeds have been specially developed for use in the new generation of 'no hospital bed' Midwife-led Birthing Centres, to provide safe, practical, comfortable and instantly accessible, rest and overnight facilities, when required, in rooms designed to appear to have no bed and look more like a smart, spacious hotel room, rather than a hospital one.

# Double Beds So Partner Can Stay Patient Overnight

When, not in use, these special counterbalanced Patient wall beds are designed to be folded away fully made-up into a wall unit. - (see also back page), but can be deployed in literally seconds, when needed. Double, rather than single beds are usually provided, so that the patient's partner can stay with them overnight — (see images 1-4).

# Fully Compliant With Multiple Additional Benefits

Our 'NEAT-A-WARD'® Patient Wallbeds are fully compliant, meeting all Infection Control, Health and Safety, Manual Handling and Flame Retardancy requirements. These beds have been specifically engineered to withstand long-term, heavy duty, daily use in Hospitals and don't just save space, but are designed to offer a range of additional benefits (to

both patient, partner and Ward staff), including safety, comfort, ease-of-handling, durability, cost and time savings and easy nursing access for staff and value for money. 'NEAT-A-WARD'® Patient Wallbeds set the standard for the sector.

# Ultra-strong Construction Guaranteed for 10 years

'NEAT-A-WARD'® Patient Wallbeds are made from hand-welded, powder-coated steel, with a unique, pre-tensioned woven steel mesh base (no wooden slats to break!). The sleeping height can be varied, depending on the length of the front and back legs specified, to allow for easier nursing access — (see also image 4). Designed for a 25 year-plus service life, subject to periodic servicing (we recommend just once every 5 years) 'NEAT-A-WARD'® Patient Wallbed frames will comfortably cope with loads of up to 200Kg (singles) and 400kgs (doubles) and are guaranteed for a full 10 years, but will last much longer.

# Easy to Clean

Our 'NEAT-A-WARD'® Patient Wallbeds and wipe-clean/waterproof 'MEDIFLEX'® Mattresses with removeable (zip-on/zip-off) covers, can be cleaned and disinfected, just like conventional Hospital beds and Mattresses, using the same (your usual) cleaning agents.



# HIGH PERFORMANCE 'NEAT-A-WARD'® PATIENT WALLBEDS



# **Multiple Safety Features**

Safety features include full counterbalancing, automatically opening and retracting front legs and, because these beds will be used by a constantly changing stream of visitors and staff, step-by-step, wipe-clean, laminated safety/operating instructions, permanently fixed to and prominently displayed on every bed.

# Choice of Sizes and Bedframe Colours and Handy Back Rest

"NEAT-A-WARD'® Patient Wallbeds are available in single, double and made-to-measure sizes, in a choice of 10 bedframe colours and, for added comfort and convenience, come complete with a back rest, which deploys and retracts automatically as the bed is lowered and raised and allows the user(s) to sit up in bed and read, or watch TV. There is a choice of a solid back rest - (images 1-2) - or a tubular steel one - (see image 4).

# Wipe-clean/waterproof 'MEDIFLEX'® Mattresses

Our high specification wipe-clean/waterproof 'MEDIFLEX'® foam mattresses are flame retardant to ignition source 5, come with zip-on/zip-off replaceable covers and can be cleaned just like a conventional hospital mattress using your usual cleaning agents.

# FREE Pre-sales Advisory Service

We can provide FREE expert pre-sales advice to Hospitals, PFI contractors and Architects, regarding the best model(s)/size(s) and sleeping heights for the project in-hand and can also supply beds for mock-up rooms, on larger projects.

# Full Spares and Aftersales Service

We provide full aftersales support, including spare parts/repairs for the bedframes, an optional installation service and optional periodic servicing (we recommend just once every 5 years) to help extend the already impressively long service lives of our Wallbeds.

# Easy and Intuitive to Use

'NEAT-A-WARD'® Patient Wallbeds come with a really comfortable wipe-clean/waterproof 'MEDIFLEX'® mattress, are designed to be

folded-away fully made-up (with all the bedding, including the pillows), are easy and intuitive to operate and, being counterbalanced and fitted with automatic leg action and can be folded away effortlessly and in literally seconds - (in an emergency, for instance). Special wipe-clean, flexible, "QUIK-FIT'®, bedding retaining straps - (see also image 4) - hold the mattress and all the bedding in place, when the bed is folded away.

### Benefits for Patients/Partners

The 'NEAT-A-WARD'® Patient Wallbeds can be ready to use in seconds (just open the doors, lower the bed and release and store the 'QUIK-FIT'® bedding retaining straps, safely out of the way), provides a high degree of comfort, is suitable for both short, or long stays and can be operated by the patient, without the need to call a member of staff, as it comes with full safety/user instructions.

### **Benefits for Staff**

There are significant manual handling and time saving benefits, as a result of not having to trundle temporary beds to and fro and make them up/ unmake them. Laundry bills are also reduced, due to the bedding staying on the bed and off the floor. A convenient sleeping/nursing height can be specified. A further advantage is the neat appearance of the rooms, during the day.

Because these beds can be folded away with all the bedding, they can be made-up by the housekeeping staff, during quieter periods between patients, then folded away ready for immediate use, if required, by the next patient.

# Real Value for Money

Subject to being properly maintained, 'NEAT-A-WARD'® Patient Wallbeds will give 25 years-plus of trouble-free service. This means that these beds do not have to be replaced for a good quarter of a century and represents a massive cost saving! So, 'NEAT-A-WARD'® Patient Wallbeds really are an investment rather than a cost.

water the transmission	<b>-</b>	The state of the s	ONS FOR TH	E 'NEAT-A	-WARD'® I	90* & 200H		
Sleeping (mattress) widths	(2ft.3ins) 700mm	(2ft.6ins) 750mm	(3ft. / 3ftXL) 900mm	(4ft.) 1200mm	(4ft.6ins) 1350mm	(5ft.) 1500mm	(6ft.) 1800mm	
Overall wallbed widths	705mm	780mm	935mm	1240mm	1390mm	1530mm	1870mm	
The 'NEAT-A-WARD'®	Patient Wallbed i	is available as two le	ngths: 1900mm, the	'NEAT-A-WARD	'® 190* and 2000	mm, the 'NEAT-A-	WARD'® 200#1.	
Sleeping (mattress) length*	1900mm	0mm *Standard bed dimensions detailed on side elevation drawing below.						
Sleeping (mattress) length #1.	2000mm	## Add 100mm (4ins.) to overall bed length AND height as detailed on side elevation drawing below.						
Housing Dimensions	Dimensions  For a suitable housing (with independent doors) all 500mm (20ins.) for carcass depth and ideally add 150mm (6ins.) to overall bed width for internal carcass width and 100-150mm (4-6ins.) to overall bed height # for internal carcas height.							
sleeping h of the froi folded aw housing, tl height of t	eight can be s nt and back leg ay). If it is inter ne back legs ca	beds individually specified by incre gs (this will also a nded to mount to an be made shor ee also side eleva ack page).	asing/decreasing affect the overall he bed on a plint ter than the fron	the length height when h inside the t legs by the	450		1900 · · · · · · · · · · · · · · · · · ·	,




# HIGH PERFORMANCE 'NEAT-A-WARD'® PATIENT WALLBEDS



# Wall Units to House the 'NEAT-A-WARD'® Patient Wallbeds

\*The wall units are usually architect designed, for a particular project see typical drawings below - (image 5, copyright the Frank Whittle Partnership Group of units designed for the midwife-led birthing units at Preston Royal Infirmary) - and installed by the contractor fitting out the Ward.

In addition to housing the wall bed, they provide storage space and house various medical equipment and services. To deploy the bed, first open the two bifold doors at 90 degrees to the housing and push them into the 'pockets' on either side of the bed then lower the bed. The front legs and back rest will deploy automatically.









# Housings for Smaller Projects

For smaller projects and or where architect designed housings are not appropriate we are able to supply an alternative model complete with integral housing and optional matching furniture in the form of our 'WISKAWAY'® Wallbeds scan the QR code above for further details - (see also images 6 & 7) - or call our Helpline 020 7439 1066. been able to design out any minor problems, over the years, why our his

been able to design out any minor problems, over the years, why our high performance products are a 'perfect fit' for their market and why none of our competitors' products can even get close.

# Why our products are unbeatable

Unlike most other wall bed manufacturers, who only see their products when they leave the factory, looking 'new and shiny', we maintain close contact with our customers, so we know exactly what our Visitor beds look like, after 5,10, 20 and even 30 years of active service. In conjunction with our policy of continuous product improvement, this is how we have been able to design out any minor problems, over the years, why our high performance products are a 'perfect fit' for their market and why none of our competitors' products can even get close.